Revestimiento SF (superfine)

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UNIVERSAL PHOSPHATED INVESTMENT FOR ALL KIND OF ALLOYS FOR CROWNS AND BRIDGES.
Type 1 – Class 1

IDEAL WORKING TEMPERATURE: Powder and liquid between 20 and 23 °C

POWDER-LIQUID RATIO: 100 g / 24 ml

1. MIXING AND POURING

IMPORTANT: **Revestimiento SF (superfine)** Must not come into contact with gypsum or investment material containing it. Use separate instruments. If a surface tension reducing agent is used, it must be thoroughly dried in a gentle stream of air (not compressed) before pouring the investment material. The powder must be added to the liquid (pure or diluted with distilled water depending on the alloy to be used) and vigorously mixed with a spatula for 20 seconds. It is then mixed under vacuum for 45 seconds and finally vibrated for 20 seconds.

If the mixture is only made by hand, mix with the spatula for 60 seconds approximately with the cup supported on the vibrator.

If rigid rings are used, its inner wall must be coated with a fire-proof tape with one turn for the rings N° 1X y 3 X and with 2 turns for the rings 6X and 9X.

2. WORKING TIME

Approximately 6 minutes

3. SETTING TIME BEFORE PREHEATING

From 45 to 60 minutes depending on the ring size.

Revestimiento SF (superfine) can be used without a metal ring up to and including size n° 6X (free expansion).

4. PREHEATING

Put the ring in the cold furnace and slowly raise the temperature (7 to 9 °C per minute). Keep the temperature constant for 30 to 60 minutes when the furnace reaches 280°C and also on reaching the final temperature (which will depend on the metal to be cast). For this type if phosphated investment material, it is good practice to keep the temperature constant for 45 to 60 minutes when it reaches 580 °C.

5. MAXIMUM TEMPERATURE: 1050 °C

6. CASTING

The time elapsing from withdrawing the ring from the furnace to casting must be as short as possible. After casting, it must be allowed to cool to ambient temperature before mould release.

7. DIMENSIONAL CONTROL

Alloys	Protheses	Concentration of the liquid
Non-precious	Large Bridges	90 - 100%
	Crowns and small bridges	70 - 90%
	Studs, inlays and conical crowns	55 - 65%
With high PD content	Ceramic-metallics	65 - 75%
Precious and Semi-precious	Crowns and bridges	60 - 65%
	Studs, inlays and conical crowns	50 - 55%

Expansion may be reduced by diluting the special liquid. This must be done using distilled water.

8. RECOMMENDATIONS

- Store the **Revestimiento SF (superfine)** powder in a cool dry place.
- Do not store the special liquid at temperatures lower than + 6 °C
- Bear in mind that the initial setting time is shortened when the temperature is more than 23 °C


9. SAFETY INSTRUCTIONS

- This material contains silicon and inhaling it may cause lung damage (silicosis or cancer).
- Avoid inhaling the dust and use a suitable face mask for dust protection.
- Do not open the furnace during the ring preheating stage, since the wax vapours may ignite in contact with air.

10. PROPERTIES

- Fluidity: 60 - 90 dPa.
- Initial setting time: 6 minutes.
- Compression resistance: 10 - 15 MPa.
- Thermal linear temperature variance: 1.4%

We guarantee the optimum quality of our products, but we cannot be held liable for the results of their subsequent use.

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