# ventura top Light

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
## INSTRUCTIONS FOR USE

### DOSAGE

Extrude the required quantity, observing the dosage scale on the mixing block. Then dose a quantity of Catalyst Gel **twice the length of the fluid**.

### MIXING

Mix energetically with the spatula, and then carefully spread the mixture out, pressing it over the mixing block to eliminate any air which may have been incorporated into the mixture. Repeat this operation until the mix is even in colour. Mixing time should ideally be about 30 seconds.

Apply the mixture to the impression tray with a spatula or a syringe for elastomers. Working time including mixing time is about 1 minute and 30 seconds.

Put into the oral cavity and wait for the mix to harden, which will take about 3 minutes and 30 seconds. The total setting time is 5 minutes.

### WARNING

In sensitive subjects, polysiloxanes may cause irritation or other allergic reactions.

### CLEANING AND DISINFECTING THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected. When using disinfectants, always follow the manufacturer's instructions.

### CASTING MODELS

The impression must be cast in between 30 minutes to 72 hours after hardening. The material is compatible with all quality commercially available chalks. We recommend the following gypsum: Type IV **Ventura Superdie Rock** and Type III **Ventura Stone**.

### WARNING / CAUTION

- Before making the corrective impression, remove all remaining traces of oral fluids by rinsing thoroughly and drying.
- The times given are for mix starting at 23°C. High temperatures or overuse of catalysts will speed up hardening, whereas low temperatures or insufficient use of catalysts will slow it down.
- Close all hardening substances properly after use.
- Prevent catalyst from coming into contact with skin and eyes; wash well with soap and water in the event of accidental skin contact. In the event of eye contact, rinse thoroughly with running water and seek medical attention. If substance is ingested, seek medical attention immediately.
- Do not stain garments.

### GALVANISATION

The impressions may be galvanised with copper or silver.

### CONSERVATION

Guaranteed for 36 months at 5° to 27°C (41° – 80°F).

### FOR DENTAL USE ONLY

Instructions for use issued in April 2016 (2016-04)

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