# Revestimiento HS (high speed)

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## RAPID PHOSPHATED INVESTMENT FOR ALL KIND OF ALLOYS FOR CROWNS AND BRIDGES.

### Type 1 – Class 2

Revestimiento HS High Speed is a universal carbon-free investment material, suitable for all the range of alloys used in fixed prostheses

#### 1. MIXING

Place the liquid in the cup first and then slowly add the powder so that it is uniformly wetted. If it is mixed manually, it is recommended to do so over a vibrator.

#### 2. MIXING RATIO AND TIME

Powder/liquid ratio: 100 g/ 24-26 ml.  
-Mix manually for 30 seconds.  
-Mix 1 minute under vacuum.  
-Keep under vacuum for a further 30 seconds.

#### 3. WORKING TIME

Approximately 5 minutes when the ambient temperature, the powder and the liquid are at a temperature between 19 and 23 °C. With the liquid hotter the working time is less and prolonged when it is cold.

#### 4. CASTING

Apply a surface tension reducer to the wax and then pour slowly the investment material on one side of the ring wall.

If steel rings are used, their inner wall must be covered with a fire-proof tape (asbestos substitute) in the following way:

- Rings No. 1 or 3: one turn of the tape
- Rings No. 6 or 9: two turns of the tape

When steel rings are not used, rings of plastic or other materials may be used that have their inner wall smooth. For those cases where the edges overlap, e.g., magnetic rings, vaseline should be applied.

For the large size No. 9 X, it is always necessary to work with metal rings.

#### 5. TIME BEFORE PREHEATING

Some 30 minutes after preparing the mixture, the ring can be placed in the already preheated furnace at a temperature of between 700 and 800 °C for conventional precious metal alloys, between 850 and 900 °C for precious metal alloys for ceramic and between 900 and 950 °C for non-precious metal alloys.

#### 6. PREHEATING TIMES

Depending on the ring size, the time left at the final temperature must be:  
-Ring No. 1X: 25 minutes.  
-Ring No. 6X: 45 minutes.  
-Ring No. 3X: 35 minutes.  
-Ring No. 9X: 60 minutes

#### 7. FINAL TEMPERATURES

Although the investment material can be heated to a maximum temperature of 1050 °C, the final pre-heating temperature of the ring depends on the alloy to be cast and is that indicated in the previous section 5 "PREHEATING TEMPERATURES".

#### 8. DIMENSIONAL CONTROL

Although it depends on the alloy type to be cast, there follow some reference values for diluting the special liquid with distilled water:

Special Liquid	Works
50%	Onlays, inlays y conical crowns.
65%	Crowns and bridges in conventional precious and semi-precious alloys.
75%	High palladium content alloys in the metal-ceramic technique.
100%	Non-precious metal alloys.

#### 9. MOULD RELEASE

Let the ring cool to ambient temperature.

#### 10. RECOMMENDATIONS

- Store the HS High Speed powder in a cool dry place.
- Do not store the special liquid at temperatures below +6 °C.
- Clean the mixing pumps and dosers well, avoiding their coming into contact with gypsum.


#### 11. SAFETY INSTRUCTIONS

- This material contains silicon and inhaling it may cause lung damage (silicosis or cancer)
- Avoid inhaling the dust and use a suitable face mask for dust protection.
- Do not open the furnace in the ring preheating phase, as the wax vapours may ignite in contact with the air.

#### 12. PROPERTIES

- Fluidity: 45 - 65 dPa.
- Initial setting time: 6 minutes.
- Compression resistance: 10 - 15 MPa.
- Linear thermal dimensional variance: 1.2%

#### 13. INFORMATION OF INTEREST FOR WORKING WITH RAPID INVESTMENT

From filling the ring with investment material, the flowing times must elapse (depending on ring size) before placing them in the oven already preheated to 850 - 900 °C

- 1 X: 20 minutes.
- 3 X: 20-23 minutes.
- 6 X: 23-28/30 minutes.

It is not recommended to work with the rapid version when more than 30 minutes have elapsed from filling the ring with investment material. In this case, the traditional method must be used.

When a silicone ring is used, great care must be taken during the setting time, avoiding pressing or pushing the ring during this time with the hands or pincers.

#### 14. ADDITIONAL INSTRUCTIONS FOR WORKING WITH INJECTED CERAMICS

Expansion liquid concentration

- For inlays - onlays (50 - 55%)
- For veneers: 60%
- For simple crowns: 60 - 65%
- For small bridges (maximum 3 elements): máx. 70%

Rings for injected ceramic must remain in the furnace at least the following times:

- Small rings: 1 X: 50 minutes
- Medium rings: 3 X: 60 minutes

For each additional ring, the time in the furnace must be increased by ten minutes.

For mould release of the injected ceramic, allow to cool to ambient temperature and observe the ceramic manufacturer's instructions for use.

Revestimiento HS High Speed material can be used with the following ceramics: Empress I\* (Ivoclar)\*, Empress II\* (Ivoclar)\*, Authentic\* (Ceramay)\*, Generic\* (Pentron)\*, Finesse\* (Dentsply)\*. Empress I\*, Empress II\*, Ivoclar\*, Authentic\*, Ceramay\*, Generic\*, Pentron\*, Finesse\*, Dentsply\* are not registered trade marks of Madespa, S.A.

*We guarantee the optimum quality of our products, but cannot be held liable for their subsequent production.*

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