ventura gingival

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DESCRIPCIÓN

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EN INSTRUCTIONS FOR USE:

DESCRIPTION

Ventura gingival is a medium viscosity vinyl polysiloxane materials suitable for the reproduction of gingival masks on master models. They are precise silicones which keep their accuracy for days due to their outstanding dimensional stability. The highest aesthetic results can be obtained because of their natural gingival colour. For extra-oral dental use only.

TECHNICAL DATA

Quick and accurate gingival reproduction on master models.

WORKING AND SETTING TIMES

Mixing ratio	1:1
Mixing time	auto-mix system
Mixing tip	Rigid: Green Elastic: Yellow
Total working time	1' 00"
Setting time	5' 00"

ATTENTION!

Application is at 23° C ± 2° C / 73 ° F ± 4 ° F, 50 ± 5% relative humidity. Higher temperature reduces working and setting times and lower temperature prolongs working and setting times.

APPLICATION

Mounting the cartridge

Insert cartridge onto dispensing gun following the manufacturer's instructions and remove cap. Extrude impression material until it exits both chambers at the same rate and then mount the appropriate mixing tip.

Application

Apply **ventura gingival** according to the desired technique (direct or indirect).

Extrude material keeping the end of the mixing tip immersed in the material to avoid the formation of air bubbles. Carefully remove the set material after the recommended time. After use keep mixing tip on cartridge as it serves as a cap. Mixing tips are for single use only.

TECHNICAL DATA

	ELASTIC	RIGID
Recovery from	>99.8	>99.8
Linear dimensional	<0.3	<0.3
Hardness Shore A	47	75
Tear strength (N/mm)	>6.5	>6.5

SCANNING

Ventura gingival elastic can be scanned with the help of a scanning spray. **Ventura gingival rigid** can be scanned with most optical and laser scanners.

PRECAUTIONS

Contact with latex gloves may impair the setting reaction. Vinyl gloves are recommended.

Do not combine with condensation curing silicones.

Avoid contact of materials with eyes. In case of contact rinse eyes with plenty of water and consult a physician immediately.

Avoid spillage of materials on clothing as they cannot be removed chemically.

STORAGE

Store material between 10-25°C (50-77°F) Store impression at room temperature (max. 25°C/77°F).

LOT NUMBER AND EXPIRY DATE

The above are printed on each package and should be quoted in all correspondence which requires identification of the product.

FOR PROFESSIONAL DENTAL LABORATORY USE ONLY

These materials have been developed for use exclusively by laboratory technicians and must be handled according to the instructions of use. The recommendations included in the instructions correspond to the actual state of knowledge regarding dental techniques. The manufacturer cannot be held responsible for damages caused by other uses or incorrect handling.

SEPARATOR FOR VENTURA GINGIVAL

Ventura separator for **ventura gingival** is a solution for the prevention of adhesion between compatible laboratory silicones in impressions or silicone moulds. **Ventura separator** can be applied to reproduce gingiva on models in association with **ventura gingival**.

INDICATIONS

- Separator for C-Silicones (polysiloxanes)
- Separator for A-Silicones (vinylpolysiloxanes)

APPLICATION

This product can be applied either directly onto the dry impression (direct technique) or onto silicone moulds (indirect technique).

Make sure that the surfaces to be isolated are thoroughly cleaned and dried. Spread Separator using the brush.

Wait for Separator to dry completely, until the surface becomes opaque. Drying time can be shortened by using a gentle blow of air.

PRECAUTIONS

- Shake well before use
- Close cap immediately after use

Instructions for use issued in February 2018 (2018-02).

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