

# NEW ventura top 2.0

*VENTURA TOP condensation silicone is now VENTURA TOP 2.0!*



## VENTURA QUALITY



*The formula is practically the same, and both its use for impression taking and handling is the same as the previous version.*

## NEW COLOURS

*The main change is in the color of the product:*

**COD. 3388 Ventura Top 2.0 Catalyst:** now **deep blue**.

*The previous version was red.*

**COD. 3387 Ventura Top 2.0 Light:** now **light pink**.

*The previous version was blue.*

**COD 3386 Ventura Top 2.0 Putty:** now **light blue**.

*The previous version was yellow.*



## IDEAL USE

*It is possible to combine the previous version with the new Top 2.0 since they are completely compatible, both in the mixture with the catalyst and in the impression. It does not alter the result.*

*It should be noted that the new Top 2.0 has a longer shelf life than the previous version, from 3 to 4 years.*

